

Products	My Own Skin	Avita ReCell	Polarity Skin TE	Epicel
Q Code	Q4226	None	Q4200	Q4100
Applications	Suture/ Staple/Glue	Spray	Slurry	Staple
Autologous	YES	YES	YES	YES
100% Autologous	YES	NO	NO	NO
Tensile Strength	3-5 times an autograft	N/A	N/A	N/A
Sizes	25cm/50 cm/ 100cm sheets	N/A	N/A	50 cm
Biopsy Size	1 cm	1cm to 24cm	3 cm to 24 cm	6cm x 2 cm
Days to Make	7 days	2 hours	2-3 days	19 days
Product Thickness	2mm	N/A	N/A	2-8 cells
Donate blood, Plasma after	YES	YES	YES	NO
Applications	Normally ONE	MULTIPLE	Normally ONE	Normally ONE